

PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT  
FOR FUTURE USE IN CASE YOU HAVE A PROBLEM.

**BAUSCH  
& LOMB**

**ReNu**®

#### Rewetting Drops

**Bausch & Lomb® ReNu** Rewetting Drops is a sterile, preserved solution recommended to effectively lubricate, rewet and help prevent film deposit buildup on soft contact lenses (including frequent replacement and disposable lenses) for clearer, more comfortable lens wear.

**CONTENTS:** **Bausch & Lomb ReNu** Rewetting Drops is a sterile, isotonic solution that contains boric acid, poloxamine, sodium borate, and sodium chloride; preserved with sorbic acid (0.15%) and edetate disodium (0.1%).

**ACTIONS:** **Bausch & Lomb ReNu** Rewetting Drops rewets lenses and helps prevent film deposit buildup. Ideal for rewetting extended wear, frequent replacement or disposable lenses, especially after awakening. Also, helps remove particulate matter.

**INDICATIONS (Uses):** Use to lubricate and rewet your lenses during wear. If you wear extended wear, frequent replacement or disposable lenses, this solution can be used to rewet your lenses after sleep.

**CONTRAINDICATIONS (Reasons not to use):** If you are allergic to any ingredients in **Bausch & Lomb ReNu** Rewetting Drops, do not use this product.

**WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.** It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping.

Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight.

Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner.

It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

- To avoid contamination, do not touch tip of container to any surface.
- Replace cap after using.

#### PRECAUTIONS:

- Always clean, disinfect, rinse and soak your lenses each time they are removed.
- Keep the bottle tightly closed when not in use.
- Store solution at room temperature.
- Use before the expiration date marked on the bottle and the carton.

#### ADVERSE REACTIONS (Problems and what to do):

The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

#### • IMMEDIATELY REMOVE YOUR LENSES.

- If the discomfort or problem stops, then look closely at the lens for dirt, an eyelash, or other foreign body on it, or damage such as a rough area, a tear, coatings, or heavy deposits.
- If the lens is in any way damaged, **DO NOT** put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens is a non-disposable or frequent replacement lens and has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it. If the problem continues, **IMMEDIATELY** remove the lens and consult your eye care practitioner.
- If the lens is a disposable lens and has dirt, an eyelash, or other foreign body on it, or the problem stops, discard the lens and replace it with a new disposable lens. If the problem continues after removal or upon insertion of a new lens, **IMMEDIATELY** remove the lens and consult your eye care practitioner.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

#### GOOD LENS CARE PRACTICES:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them if you are wearing non-disposable or frequent replacement lenses.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse lens case with fresh rinsing solution, and allow to air dry.
- When wearing disposable lenses, each time you remove a lens, discard it and replace it with a new lens.

#### DIRECTIONS:

- During wear, put two or three drops directly onto each lens whenever needed. Blink several times and wait a few moments.
- If the lens still does not feel comfortable, add another drop.
- If discomfort persists, **IMMEDIATELY** remove the lens and contact your eye care practitioner.

**HOW SUPPLIED:** **Bausch & Lomb ReNu** Rewetting Drops is available for purchase in 0.5 Fl. Oz. (15ml) plastic bottles. Bottles and cartons are marked with a lot number and expiration date.

**Bausch & Lomb** and **ReNu** are registered trademarks of Bausch & Lomb Incorporated.

© Bausch & Lomb Incorporated

Manufactured by:  
Bausch & Lomb, Rochester, NY 14692-0450. Made in USA.